

The New Standard in Patient Comfort

eyecia II



eyecia II


More comfort. Better care.

The EYECIA II is the quietest, most comfortable chair we've ever designed.

So even the most fearful patients can sit back and relax - and doctors can provide superior care.



EYECIA II is ergonomically designed to maximize patient comfort and relaxation. For doctors, we've simplified the controls, making key functions easier to find and use.



Designed with patients and doctors in mind.

※Image is optional dental chair.

Serenity



Smooth Start & Stop Technology

The chair starts and stops so smoothly and slowly, patients hardly feel the chair moving at all.

Silent Movement System

No jolts or mechanical sounds. The chair's motor is whisper-quiet, so patients don't hear a thing.



Premium Comfort


Our optional Takumi seat is like a custom-made lounge chair with deep, supportive cushioning.

Color options

Standard



Soft, silky upholstery that's easy to clean and makes patients feel like they're lying on a soft blanket.



Takumi seat

optional



Exceptionally soft upholstery, crafted by hand. This optional chair is like a luxury sofa. It's so comfortable, your patients will remain at ease for the entire treatment.



Hygienic

Removable Spittoon Bowl

Completely removable glass bowl is easy to clean, giving your patients peace of mind.



Automatic Spittoon Flush & Cup Filler

Convenient auto-control system for spittoon flush and cup filler - equipped as a standard.



Removable Table Tray

Removable waterproof table tray fits perfectly into the table to keep liquids out.



Touchelss Sensor Operating Light

Natural brightness and a shadow-free light field offer optimal vision. Staff can turn the light on/off hygienically without touching it.

Simplicity



*Intuitive control panel
simplifies operation*

A smarter layout of controls enables doctors to concentrate more on the procedure and on conversations with patients.



Efficiency

The most frequently used functions are all located on the front control panel. Individual instrument controls are located on the side for enhanced efficiency.



Easy Handling

The handles on both sides allow for easy repositioning of the table.



User-Friendly

The instrument holder is tilted at a 20-degree angle, improving accessibility and protecting doctors' hands from other instruments.

Utility



Tools for Every Need

A wide selection of optional devices and instruments are available to meet a variety of clinical needs.

YOSHIDA INSTRUMENTS

※ WITH OPTICAL LIGHT



V type syringe



Multiway syringe



Ultrasonic scaler

AIR TURBINE



Royalroad SDC1



Royalroad SDC1mini



Super load 640BP-LUX



Super load 640BQP



High torque 635BP-LUX



High torque 635BQP

MICROMOTOR HANDPIECE



S740 MI



S740 SI

QUICK JOINT




LX503(LED)



NO.3

MICROMOTOR



40C

DISPLAY ARM



Display arm & monitor

UNIT




Water bottle

FILM VIEWER



Panoramic film viewer

FOOT CONTROLLER






Step-in type



Sliding foot controller

Dimensions



CONTACT

Copyright © 2014 THE YOSHIDA DENTAL MFG. CO., LTD.

All right reserved. No part of this publication may be reproduced in any form without the written permission of The Yoshida Dental Mfg. Co., Ltd.



MD 554641/ISO 13485



FM 554640/ISO 9001

YOSHIDA

THE YOSHIDA DENTAL MFG. CO., LTD.

Address 1-3-6 Kotobashi, Sumida-ku Tokyo, Japan ZIP 130-8516

TEL +81-3-3631-2165 FAX +81-3-3631-2685 (International Business Div.)